

Genome Sequence of *Bacillus alveayuensis* Strain 24KAM51, a Halotolerant Thermophile Isolated from a Hydrothermal Vent

Sevasti Filippidou,^a Tina Wunderlin,^a Thomas Junier,^b Nicole Jeanneret,^a  Shannon Johnson,^c Kim McMurry,^c Cheryl D. Gleasner,^c Chien-Chi Lo,^c Po-E Li,^c Momchilo Vuysich,^c Patrick S. Chain,^c Pilar Junier^a

Laboratory of Microbiology, Neuchâtel, Switzerland^a; Vital-IT Group, Swiss Institute of Bioinformatics, Lausanne, Switzerland^b; Los Alamos National Laboratory, Los Alamos, New Mexico, USA^c

***Bacillus alveayuensis* strain 24KAM51 was isolated from a marine hydrothermal vent in Milos, Greece. Its genome depicts interesting features of halotolerance and resistance to heavy metals.**

Received 17 July 2015 Accepted 21 July 2015 Published 27 August 2015

Citation Filippidou S, Wunderlin T, Junier T, Jeanneret N, Johnson S, McMurry K, Gleasner CD, Lo C-C, Li P-E, Vuysich M, Chain PS, Junier P. 2015. Genome sequence of *Bacillus alveayuensis* strain 24KAM51, a halotolerant thermophile isolated from a hydrothermal vent. *Genome Announc* 3(4):e00982-15. doi:10.1128/genomeA.00982-15.

Copyright © 2015 Filippidou et al. This is an open-access article distributed under the terms of the [Creative Commons Attribution 3.0 Unported license](https://creativecommons.org/licenses/by/3.0/).

Address correspondence to Pilar Junier, pilarjunier@unine.ch.

Bacillus alveayuensis is a thermophilic endospore-forming bacterium. The type strain MT1 was isolated from deep-sea sediment (1). However, to date, no genome is available for this species. Strain 24KAM51 was isolated from a hydrothermal vent on the coastal line of Alykes beach in Milos, Greece (36°42'353" N, 24°28'197" E, depth 1.5 m). It is able to grow at a temperature range of 30 to 80°C, with an optimum at 60°C. It can also tolerate acidic and alkaline conditions (pH growth range from 3 to 10), with an optimum at 7. Finally, it grows with up to 13% (wt/vol) NaCl. Based on 16S rRNA sequence identity, strain 24KAM51 is closely related to *B. alveayuensis* type strain TM1 (99% identity). A series of physiological tests also showed similarity to TM1. The genome of 24KAM51 was sequenced and annotated in order to contribute to a better understanding of the thermophilic and halotolerant lifestyle of extremophilic bacilli. To date, this is the only genome of *B. alveayuensis* publicly available.

Genomic DNA was extracted from an overnight culture using the QIAamp DNA minikit (Qiagen GmbH, Germany). For the draft genome of *B. alveayuensis* strain 24KAM51, an Illumina (2) short-insert (300 ± 70 bp) library was constructed and sequenced generating 43,458,058 reads, totaling 4.39 Mbp. The Illumina draft data were assembled with Velvet, version 1.2.08 (3). The estimated size of the genome is 6.7 Mbp, which provides 409× coverage of the genome. Genome annotation was performed using an Ergatis-based (4) workflow with minor manual curation and visualized with the Artemis genome browser and annotation tool (5). The G+C content is 38.1%. The complete genome sequence contained 6,597 genes, 10 rRNAs (5S, 16S, and 23S), 88 tRNAs, and 2 noncoding RNAs (ncRNAs) predicted. The genome of 24KAM51 was annotated, and its proteome revealed the presence of 198 genes related to the sporulation and germination pathways and 10 genes related to dipicolinic acid synthesis. Fifty-three proteins related to flagellar motility were found in the proteome of *B. alveayuensis* strain 24KAM51.

Many proteins are related to halotolerance in bacilli (6), some

of which are present in multiple copies in the 24KAM51 proteome. More precisely, genes for Na⁺/H⁺ antiporters, proline/Na⁺ symporters, and glycine/betaine ABC transporters were found. Although isolated from a natural hot spring, its genome contains genes related to copper (copper-binding proteins and multicopper oxidase), manganese (manganese transporters and permeases), cadmium (cadmium transporter), zinc (zinc metalloproteases, proteases, and transporters), and arsenic (arsenic resistance protein, ArsB) resistance. Moreover, it possesses genes encoding the NarH and NarZ proteins (nitrate reduction), as well as DsrE (sulfur reduction).

Nucleotide sequence accession numbers. This whole-genome shotgun project has been deposited at DDBJ/EMBL/GenBank under the accession no. [JYCE00000000](https://www.ncbi.nlm.nih.gov/nuccore/JYCE00000000). The version described here is JYCE00000000.1.

ACKNOWLEDGMENTS

This work was financially supported by the Swiss National Science Foundation project 31003A_152972 and by Fondation Pierre Mercier pour la Science.

The genetic information downloaded from GenBank is considered to be part of the genetic patrimony of Greece, the country from which the sample was obtained. Users of this information agree to: 1) acknowledge Greece as the country of origin in any country where the genetic information is presented and 2) contact the CBD focal point and the ABS focal point identified at the CBD website <http://www.cbd.int/information/nfp.shtml> if they intend to use the genetic information for commercial purposes.

REFERENCES

- Bae SS, Lee J-H, Kim S-J. 2005. *Bacillus alveayuensis* sp. nov., a thermophilic bacterium isolated from deep-sea sediments of the Ayu Trough. *Int J Syst Evol Microbiol* 55:1211–1215. <http://dx.doi.org/10.1099/ijs.0.63424-0>.
- Bennett S. 2004. Solexa Ltd. *Pharmacogenomics* 5:433–438. <http://dx.doi.org/10.1517/14622416.5.4.433>.
- Zerbino DR, Birney E. 2008. Velvet: algorithms for *de novo* short read

- assembly using de Bruijn graphs. *Genome Res* 18:821–829. <http://dx.doi.org/10.1101/gr.074492.107>.
4. Hemmerich C, Buechlein A, Podicheti R, Revanna KV, Dong Q. 2010. An Ergatis-based prokaryotic genome annotation Web server. *Bioinformatics Oxf Engl* 26:1122–1124. <http://dx.doi.org/10.1093/bioinformatics/btq090>.
 5. Carver TJ, Rutherford KM, Berriman M, Rajandream M-A, Barrell BG, Parkhill J. 2005. ACT: the Artemis comparison tool. *Bioinformatics* 21:3422–3423. <http://dx.doi.org/10.1093/bioinformatics/bti553>.
 6. Takami H, Takaki Y, Uchiyama I. 2002. Genome sequence of *Oceanobacillus iheyensis* isolated from the Iheya Ridge and its unexpected adaptive capabilities to extreme environments. *Nucleic Acids Res* 30:3927–3935. <http://dx.doi.org/10.1093/nar/gkf526>.